

Chemiluminescence

Hideharu Shintani*

Chuo University, School of Science, 1-13-27, Kasuga Bunkyo 112-0003 Tokyo, Japan

Keywords: Chemiluminescence; Superoxide anion radical; Leukocytes; Myeloperoxidase

Introduction

Inflammatory cells such as neutrophils, eosinophils, and macrophages produce a multitude of reactive oxygen species (ROS) when they are activated. Of the reactive oxygen species, superoxide anion radical (O_2^-), a precursor of several oxidants such as hydrogen peroxide (H_2O_2), hydroxy radical (OH^\cdot), singlet oxygen ($^1\text{O}_2$), and hypochlorite (HOCl), can be formed in, and released from, each of three cell lines. On the other hand, $^1\text{O}_2$ (singlet oxygen, metastable oxygen) and hypohalides can be generated by myeloperoxidase (MPO)-mediated reactions in the neutrophil phagosome at pH 4.5, and by the eosinophil peroxidase-mediated reaction in eosinophils at physiological pH. It is, therefore, important to detect and measure the generation of O_2^- and $^1\text{O}_2$ in small amounts in such biological systems. Because production of highly bactericidal agents such as $^1\text{O}_2$ and HOCl is a result of myeloperoxidase activity, simple and sensitive methods to measure this enzyme activity are required. Herein, we describe the measurement of O_2^- , $^1\text{O}_2$ generation, and myeloperoxidase activity, in leukocytes by chemiluminescence using cypridina luciferin analogues as chemiluminescence probes.

O_2^- Production by Stimulated Leukocytes

When granulocytes or macrophages are stimulated, they produce and release O_2^- . The O_2^- released from these cells reacts efficiently with cypridina luciferin analogue (CLA, 2-methyl-6-phenyl-3,7-dihydroimidazo[1,2-a]pyrazin-3-one. MCLA, 2-methyl-6-(p-methoxyphenyl)-3,7-dihydroimidazo[1,2-a]pyrazin-3-one) (Figure 1) to produce an excited carbonyl which emits light in the visible region; i.e., maximum emission at 380 nm for CLA and at 457 nm for MCLA [1-3].

Protocol

1. Dissolve powdered Hanks' balanced salt solution (HBSS) for tissue culture (Nissui Pharmaceutical Co. Ltd, Tokyo) in

distilled water (9.8 g/L) and adjust to pH 7.4 by use of aqueous Na_2HPO_4 solution (280 mOsm).

2. Prepare granulocytes or macrophages by standard methods.
3. Dissolve MCLA ($E_{430\text{ nm}} = 9600\text{ M/cm}$) or CLA ($E_{410\text{ nm}} = 8900\text{ M/cm}$) in distilled water at a concentration of $40\text{ }\mu\text{M}$.
4. Dilute granulocytes or macrophages with HBSS at 5×10^4 cells for MCLA or 4×10^5 cells for CLA.
5. Add MCLA or CLA ($50\text{ }\mu\text{L}$, final concentration: $1\text{ }\mu\text{M}$).
6. Preincubate at 37°C for 3 min.
7. Add opsonized zymosan (OZ, in water, 20 mg/mL , $80\text{ }\mu\text{L}$), formyl-methionyl-leucyl-phenylalanine (fMLP in Ca^{2+} - and Mg^{2+} -free HBSS/ DMSO = 1:1, v/v, $200\text{ }\mu\text{M}$, $10\text{ }\mu\text{L}$), or 4β -phorbol myristate acetate (PMA in Ca^{2+} - and Mg^{2+} -free HBSS/DMSO = 1:1, v/v, $2\text{ }\mu\text{g/mL}$, $20\text{ }\mu\text{L}$) (Total volume 2.0 mL).
8. Measure chemiluminescence.


MPO in granulocytes $^1\text{O}_2$ measurements

When peripheral leukocytes derived from bone marrow are separated the neutrophils contain abundant granule myeloperoxidase as a host enzyme. Under acidic conditions, myeloperoxidase (MPO) catalyses the production of HOBr and HOCl in the presence of Br^- and Cl^- , respectively (Reaction 1). The HOBr or HOCl can efficiently react with H_2O_2 to generate $^1\text{O}_2$ (Reaction 2) [4].

MPO




pH 4.5-5.0




MCLA reacts not only with O_2^- but also with $^1\text{O}_2$ to emit light with maximum intensity at 457 nm . The addition of $0.5\text{ }\mu\text{M}$ SOD to the system, in which both O_2^- and $^1\text{O}_2$ are generated, completely quenches the O_2^- dependent luminescence leaving only the $^1\text{O}_2$ -dependent luminescence.

Protocol

1. Mixing human neutrophils and hexadecyltrimethylammonium bromide (0.02% w/v) in potassium phosphate buffer (pH 6.0, 50 mM) for 10 min.
2. Sonicate the mixture on ice for 30 s (0.5-s bursts at 0.5-s



A5307



A5309

Figure 1: Chemical structure of CLA and MCLA.

*Corresponding author: Hideharu Shintani, Chuo University, School of Science, 1-13-27, Kasuga Bunkyo 112-0003 Tokyo, Japan, Tel: +81425922336; E-mail: shintani@mail.hinocav.ne.jp

Received May 20, 2013; Accepted June 26, 2013; Published June 28, 2013

Citation: Shintani H (2013) Chemiluminescence. Pharm Anal Acta 4: 252. doi:10.4172/2153-2435.1000252

Copyright: © 2013 Shintani H. This is an open-access article distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited.

intervals) at an output level of 30 W (Branson Sonifier 250, Branson Ultrasonics Corporation, Danbury, CT, USA).

3. Centrifuge at 40000 g for 20 min at 4°C (The supernatant is used as the sample).
4. Preincubate a standard solution containing MCLA (in water, 200 µM, 100 µL) SOD (in water, 50 µM, 20 µL), KBr (in water, 50 µM, 20 µL), desferrioxamine (in water, 2 mM, 20 µL), H₂O₂ (40 mM, 25 µL), and acetate buffer (0.2 M, 1 mL) (Total volume 1.95 mL).
5. Measure chemiluminescence until baseline is constant.
6. Inject samples rapidly (50 µL) by means of a microsyringe.
7. Measure chemiluminescence.

References

1. Nakano M, Sugioka K, Ushijima Y, Goto T (1986) Chemiluminescence probe with Cypridina luciferin analog, 2-methyl-6-phenyl-3,7-dihydroimidazo[1,2-a]pyrazin-3-one, for estimating the ability of human granulocytes to generate O₂⁻. Anal Biochem 159: 363-369.
2. Nishida A, Kimura H, Nakano M, Goto T (1989) A sensitive and specific chemiluminescence method for estimating the ability of human granulocytes and monocytes to generate O₂⁻. Clin Chim Acta 179: 177-181.
3. Sugioka K, Nakano M, Kurashige S, Akuzawa Y, Goto T (1986) A chemiluminescent probe with a Cypridina luciferin analog, 2-methyl-6-phenyl-3,7-dihydroimidazo[1,2-a]pyrazin-3-one, specific and sensitive for O₂-production in phagocytizing macrophages. FEBS Lett 197: 27-30.
4. Uehara K, Hori K, Nakano M, Koga S (1991) Highly sensitive chemiluminescence method for determining myeloperoxidase in human polymorphonuclear leukocytes. Anal Biochem 199: 191-195.