

# Nitrosothiol Detection by HPLC Coupled with Flow Reactors of $\rm Hg^{2+}$ and Griess Reagent

## Hideharu Shintani\*

Chuo University, School of Science, 1-13-27, Kasuga Bunkyo 112-0003 Tokyo, Japan

Keywords: Nitrosothiol; HPLC; S-nitrosoproteins; Quantification

# Introduction

NO-related intermediates, including NO<sup>+</sup> like species, can participate in nitro-sating additions to nucleophilic centers of biological molecules. Sulfhydryl-containing molecules such as glutathione are particularly susceptible to nitrosation and form nitrosothiol adducts (nitrosothiols; RS-NOs). It seems that these adducts in biological systems play an important role in NO-mediated signalling cascades such as the downregulation of N-methyl-D-aspartate receptor, and the regulation of transcriptional factors [1]; they might also be involved in non-adrenergic and non-cholinergic neuronal responses [2]. It is, therefore, essential to identify nitrosothiols specifically and to quantify them in biological systems.

Herein we describe a sensitive and specific HPLC method coupled with  $Hg^{2+}$  and Griess reagent for nanomolar quantification of a wide range of RS-NOs including low-molecular weight RS-NOs such as nitrosoglutathione (GS-NO) and nitrosocysteine (Cys-NO), and also S-nitrosoproteins, in particular S-nitrosoalbumin and S-nitrosohaemoglobin [3].

## **Preparation of S-nitrosoproteins**

#### Protocol

**Reagents:** Bovine serum albumin (BSA) (Nacalai Tesque, Osaka, Japan), dithiothreitol (DTT; Wako Pure Chemicals, Osaka, Japan), isoamyl nitrite (Wako), ethyllenediaminetetraacetic acid (EDTA; Dojindo Laboratories, Kumamoto, Japan), diethylenetriaminepentaacetic acid (DTPA; Dojindo), and 5, 5'-dithiobis[2-nitrobenzoic acid] (DTNB; Nacalai).

Procedure: Reduction of BSA

- 1. Add DTT (100 mM, 10  $\mu$ L) to BSA (1 mM, 1 mL) in sodium phosphate buffer (pH 7.0, 100 mM).
- 2. Incubate for 30 min at 37°C.
- 3. Apply to a Sephadex G-100 (Pharmacia, Uppsala, Sweden) column equilibrated with sodium phosphate buffer (pH 7.0, 10 mM) containing EDTA (1 mM; buffer A), and elute with buffer A.
- 4. Collect the fractions containing BSA and concentrate the fractions by ultrafiltration.

Nitrosation of the reduced BSA: 1. Add isoamyl nitrite (100 mM, 10  $\mu$ L) to reduced BSA (0.1 mM, 1 mL) in sodium phosphate buffer (pH 7.8, 100 mM) containing DTPA (0.5 mM).

- 2. Incubate for 30 min at 37°C.
- 3. Apply to a Sephadex G-25 (Pharmacia) column equilibrated with buffer A, and elute with buffer A.
- 4. Collect the fractions containing S-NO-BSA and concentrate the fractions by ultrafiltration.

## **Quantification of Nitrosothiols**

## Protocol

**Reagents:** Nitrosoproteins (e.g., S-NO-BSA) prepared above, GS-NO (Dojindo), Cys-NO (Dojindo), HgCl<sub>2</sub>, and Griess reagent (naphthylethylenediamine; sulfanilamide (Wako)).

**Principle:** RS-NOs of different molecular size separated by HPLC are converted to NO<sup>+</sup> or NO<sub>2</sub><sup>-</sup>, in a flow reactor system connected serially to the HPLC (Figure 1), by rapid and quantitative metal-catalysed reaction with Hg<sup>2-</sup>. As reported by Saville, Hg<sup>2+</sup> decomposes the nitrosothiols stoichiometrically to nitrite according to the equations [4,5]:

 $RS-NO + Hg^{2+} \Leftrightarrow [RS-(Hg)NOj^{2+} [4]]$ 

 $[RS-(Hg)NOj^{2+} + H_2O \rightarrow RSHg^+ + NO_2^{-} + 2H^+ [5]]$ 

Peak detection is based on the colorimetric assay of nitrite using Griess reagent, which forms a diazo dye having strong absorbance at 540 nm under acidic conditions; the extinction coefficient is 53000 M/ cm [6].

#### HPLC conditions

The injection volume is 150 µL.

#### Column

1. Separation: a C1B · reversed ≈ hase column (4.6×250 mm;



\*Corresponding author: Hideharu Shintani, Chuo University, School of Science, 1-13-27, Kasuga Bunkyo 112-0003 Tokyo, Japan, Tel: +81425922336; E-mail: shintani@mail.hinocatv.ne.jp

Received May 21, 2013; Accepted June 24, 2013; Published June 28, 2013

Citation: Shintani H (2013) Nitrosothiol Detection by HPLC Coupled with Flow Reactors of  $Hg^{2*}$  and Griess Reagent. Pharm Anal Acta 4: 250. doi:10.4172/2153-2435.1000250

**Copyright:** © 2013 Shintani H. This is an open-access article distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited.

TSKgel OOS.80Ts; Tosoh, Tokyo) for low molecular weight RS-NOs; a gel filtration column (8×300 mm; Oiol-120; YMC, Kyoto) for the nitrosoproteins.

 Deproteinating column: small columns (3×10 mm) packed with C1B-based resin are placed just before the separation column and just after Hg<sup>2-</sup> flow reactor coil, in the reversedphase and gel filtration systems, respectively.

## Mobile phase

- 1. Pump 1 for HPLC: sodium acetate buffer (pH 5.5, 10 mM) containing DTPA (0.5 mM) and methanol (0-7%), for low molecular weight RS-NOs (reversed-phase HPLC); sodium acetate buffer (pH 5.5, 10 mM) containing DTPA (0.5 mM) and NaCI (150 mM) for the nitrosoprotein (gel filtration HPLC).
- 2. Pump 2 for Hg<sup>2+</sup>-flow reactor: HgCl<sub>2</sub> (1.75 mM) in sodium acetate buffer (pH 5.5, 10 mM).
- 3. Pump3forGriessreagentflowreactor:naphthylethylenediamine (0.1%) in  $H_2O$ ; sulfanilamide (1.0%) + phosphoric acid (2.0%) in  $H_2O$ .

Flow rate: 0.55 mL/min (pump 1), 0.2 mL/min (pump 2), 0.24 mL/min (pump 3)

Detector: visible detector (Eicom, Kyoto) operated at 540 nm.

Typical elution profiles of GS-NO, Cys-NO and S-NO-BSA are shown in Figure 2.

Figure 3 illustrates the correlation between the peak area of GS-NO calculated by the integrator and different concentrations examined by use of the HPLC f10w reactor. The detection limit for GS-NO was found to be >3 nM.



Figure 2: Elution profiles of different RS-NOs for the  ${\rm HPLC/Hg^2}$  flow reactor system.

A: Low molecular weight RS-NOs (reversed-phase HPLC). B: High molecular weight RS-NOs (gel filtration HPLC)



Page 2 of 2

Figure 2: Elution profiles of different RS-NOs for the HPLC/Hg<sup>2</sup> flow reactor system.

A: Low molecular weight RS-NOs (reversed-phase HPLC). B: High molecular weight RS-NOs (gel filtration HPLC)

## References

- Hausladen A, Privalle CT, Keng T, DeAngelo J, Stamler JS (1996) Nitrosative stress: activation of the transcription factor OxyR. Cell 86: 719-729.
- Rand MJ, Li CG (1995) Discrimination by the NO-trapping agent, carboxy-PTIO, between NO and the nitrergic transmitter but not between NO and EDRF. Br J Pharmacol 116: 1906-1910.
- Akaike T, Inoue K, Okamoto T, Nishino H, Otagiri M, et al. (1997) Nanomolar quantification and identification of various nitrosothiols by high performance liquid chromatography coupled with flow reactors of metals and Griess reagent. J Biochem 122: 459-466.
- Saville B (1958) A scheme for the colorimetric determination of microgram amounts of thiols. Analyst 83: 670-672.
- Williams DL (1996) S-nitrosothiols and role of metal ions in decomposition to nitric oxide. Methods Enzymol 268: 299-308.
- Wishnok JS, Glogowski JA, Tannenbaum SR (1996) Quantitation of nitrate, nitrite, and nitrosating agents. Methods Enzymol 268: 130-141.