

Docking Molecular Analysis of Potential Drug Paritaprevir against *Mycobacterium tuberculosis* (Mtb)

Ivan Vito Ferrari^{1*}, Paolo Patrizio²

¹Department of Industrial Engineering, University of Rome Tor Vergata, Rome, Italy; ²Department of Industrial Engineering, University of Padua, Febbraio, Italy

ABSTRACT

Mycobacterium tuberculosis (Mtb) is the causative agent of tuberculosis, which kills 1.8 million annually. This is an infectious disease generally affects the lungs, but can also affect other parts of the body. Mtb RNA polymerase (RNAP) is the target of the first-line antituberculosis drug Rifampin (Rif). We report first time a Potential Drug Paritaprevir against with severe infectious disease, by *in silico* approach, using Auto dock Vina and Auto dock 4 (or MGL Tool), estimated with Pyrx and AM Dock Software, calculating three different important parameters: Binding Affinity (kcal/mol), estimated Ki (in nM units) and Ligand Efficiency (L.E. in kcal/mol). After a selective analysis of over 1000 drugs, processed with Pyrx (a Virtual Screening software for Computational Drug Discovery) in the Ligand Binding site pocket of the protein (ID PDB 5UHB chain C, DNA-directed RNA polymerase subunit beta), we noticed high values of these 3 parameters mentioned above of Paritaprevir, concluding that it could be an excellent candidate drug for this type of infection. Indeed, from the results of Auto dock Vina and Auto dock 4 (or Auto dock 4.2), implemented with Lamarckian genetic algorithm, LGA, through AM Dock Software, This oral drug, approved by FDA in 2014, both by Auto dock Vina and Auto dock Vina 4 has excellent Binding affinity value, ca. -10.00 kcal/mol, a Ki value 40 nM and Ligand efficiency ca-0.15 kcal/mol. These results are comparable to the drug crystallized in the above-mentioned protein, currently used against TBC.

Keywords: *Mycobacterium tuberculosis* (Mtb); tuberculosis (TBC); Auto dock 4; Paritaprevir

INTRODUCTION

Mycobacterium tuberculosis (Mtb) is the causative agent of tuberculosis, which kills 1.8 million annually. Mtb RNA polymerase (RNAP) is the target of the first-line anti tuberculosis drug Rifampin (Rif) [1]. Most infections show no symptoms, in which case it is known as latent tuberculosis. The classic symptoms of active TB are a chronic cough with blood-containing mucus, fever, night sweats, and weight loss [2]. A total of 1.4 million people died from TB in 2019 (including 208 000 people with HIV). Worldwide, TBC is one of the top 10 causes of death and the leading cause from a single infectious agent (above HIV/AIDS) [2].

In this short communication, we investigated about 1000 drugs, through *in silico* Docking approach, downloaded from Pubchem

Database (<https://pubchem.ncbi.nlm.nih.gov/>). Molecular docking has become a powerful tool for drug discovery and optimization. General speaking, there are several protein-ligand docking programs currently available. This bioinformatics methodology is one of the most frequently used methods in structure-based drug design, due to its ability to predict the binding-conformation of small molecule ligands to the appropriate target (or receptor) binding site.

In this paper particular attention, we have focused on Autodock Vina, estimated with Pyrx software, a simple Virtual Screening library software for Computational Drug Discovery (<https://pyrx.sourceforge.io/>), based on prediction the binding orientation and affinity of a ligand [3] and that of Autodock 4, opened by AMDOCK software, (a versatile graphical tool for assisting molecular docking with Autodock Vina and Autodock

Correspondence to: Ivan Vito Ferrari, Department of Industrial Engineering, University of Rome Tor Vergata, Rome, Italy, Tel/Fax: +39 3932352697; E-mail: Ivan.Vito.Ferrari@uniroma2.it

Received: August 05, 2021; **Accepted:** August 19, 2021; **Published:** August 26, 2021

Citation: Ferrari IV, Patrizio P (2021) Docking Molecular Analysis of Potential Drug Paritaprevir against *Mycobacterium tuberculosis* (Mtb). Blood Disord Transfus. 12: 473

Copyright: © 2021 Ferrari IV, et al. This is an open-access article distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited.

4, characterized to be automated docking of ligand to macromolecule by Lamarckian Genetic Algorithm and Empirical Free Energy Scoring Function [4]. Both Scoring Functions were supplied by the Scripps research institute [3-4].

MATERIALS AND METHODS

Protein and ligand preparation before docking

5UHB [Crystal structure of *Mycobacterium tuberculosis* transcription initiation complex in complex with Rifampin] was prepared manually using different software, before molecular docking analysis.

The first step was downloaded from Protein Data Bank, (<https://www.rcsb.org/structure/5UHB>) and save in pdb format. The second step was removed all unnecessary docking chains. In fact, only the C chain (DNA-directed RNA polymerase subunit beta) has been maintained and re-saved in pdb format (Figure 1). Next, were the removal of ligands and water molecules crystallized using Chimera software [5]. Later, polar hydrogens and Kollmann charges were added with Mgl Tool, (or called Auto Dock tools, a software developed at the Molecular Graphics Laboratory (MGL) [6] As a last step they were added to the protein, any missing amino acids and the whole protein was minimized with the Swiss PDB Viewer Software [7]. In this way the protein was saved in pdb format and ready for docking, through Autodock Vina, estimated with a soft virtual screening library called Pyrx. On the other hand for Ligand Preparation, the first step was to separate the crystallized drug from its protein, manually add all the hydrogens and their charges (with the MGL tool software) and minimize it with MMFF94 force field [8,9], opened Pyrx software (<https://pyrx.sourceforge.io/>). This crystallized drug, Rifampin (RFP) was docked in the same binding area in order to accurately evaluate its Binding Affinity Vina Score with its protein and if there is overlapping of the drug Rifampin, between its crystallized version and in its docking version (Figure 2). In the area of the drug crystallized with the protein's C chain, more than 1000 dockings were carried out, through Autodock Vina, estimated with Pyrx software, estimating their Binding Affinity (kcal/mol). Parameters Grid Box for Docking in Ligand Binding Site Pocket for Repeatability Binding Affinity by AM Dock Software calculated with Autodock Vina and Autodock Vina: -ID PDB 5UHB Chain: Center X (=163.84); Centre Y (=162); Centre Z (=20); Dimensions (Angstrom) (Å) X, Y, Z [=17=,17, =17]; exhaustiveness =8.




Figure 1: ID PDB 5UHB Chain C (DNA-directed RNA polymerase subunit beta), comparison 3D structure cristalized and 2D plot diagram cristalized Rifampin, in ligand binding site pocket by discovery studio biovia software [9].

J Blood Disord Transfus, Vol.12 Iss.7 No:1000473




Figure 2: Comparison 3D structure and 2D-diagram, In ligand binding site pocket prepared by Discovery Studio Biovia, [9] and Chimera Software [5] in ID PDB 5UHB Chain C protein (NA-directed RNA polymerase subunit beta) complexed with crystallized Rifampin with docked Rifampin.

DISCUSSION AND RESULTS

We report first time a Potential Drug Paritaprevir, against with severe infectious disease TBC, by *in silico* approach, using Autodock Vina and Autodock 4 (or MGL Tool), estimated with Pyrx and AMDock Software, calculating three different important parameters: Binding Affinity (kcal/mol), estimated K_i (in nM units) and Ligand Efficiency (L.E. in kcal/mol). After a selective analysis of over 1000 drugs, processed with Pyrx (a Virtual Screening software for Computational Drug Discovery) in the Ligand Binding site pocket of the protein (ID PDB 5UHB chain C, DNA-directed RNA polymerase subunit beta), we noticed high values of these 3 parameters mentioned above of Paritaprevir, concluding that it could be an excellent candidate drug for this type of infection (Figure 3). Indeed, from the results of Autodock Vina and Autodock Vina 4 (or Autodock 4.2), implemented with Lamarckian genetic algorithm, LGA, trough AMDock Software, This oral drug, approved by FDA in 2014, both by Autodock Vina and Autodock 4 has excellent Binding affinity value, ca. -10.00 kcal/mol, a K_i value 40 nM and Ligand efficiency ca -0.15 kcal/mol. These results are comparable to the drug crystallized in the above-mentioned protein, currently used against TBC (Figure 3) (Table 1).




Figure 3: ID PDB 5UHB Chain C protein (DNA-directed RNA polymerase subunit beta): comparison 3D structure docked and 2D plot diagram docked Paritaprevir In ligand binding site pocket by Discovery Studio Biovia Visualizer Software [9].

Drug	Scoring function	Binding affinity (Kcal/mol)	Estimation KI Nm	Ligand efficiency (Kcal/mol)
Crystalized Rifampin	Autodock Vina	8,9	299,41	0,15
Paritaprevir	Autodock Vina	10	46,77	0,18
Crystalized Rifampin	Autodock Vina 4 (MGL tool)	9,87	58,24	0,17
Paritaprevir	Autodock Vina 4 (MGL tool)	10,08	40,86	0,18

Table 1: Comparison results of crystalized drug Rifampin and proposed drug Paritaprevir calculated Autodock Vina and Autodock 4 (or MGL Tool), estimated by AMDock Software [4].

CONCLUSION

We report first time a Potential Drug Paritaprevir against with severe infectious disease, by *in silico* approach, using Autodock Vina and Autodock 4 (or MGL Tool), estimated with Pyrx and AMDock Software, calculating three different important parameters: Binding Affinity (kcal/mol), estimated Ki (in nM units) and Ligand Efficiency (L.E. in kcal/mol). After a selective analysis of over 1000 drugs, processed with Pyrx (a Virtual Screening software for Computational Drug Discovery) in the Ligand Binding site pocket of the protein (ID PDB 5UHB chain C, DNA-directed RNA polymerase subunit beta), we noticed high values of these 3 parameters mentioned above of Paritaprevir, concluding that it could be an excellent candidate drug for this type of infection. Indeed, from the results of Autodock Vina and Autodock 4 (or Autodock 4.2), implemented with Lamarckian genetic algorithm, LGA, through AMDock Software, This oral drug, approved by FDA in 2014, both by Autodock Vina and Autodock 4 has excellent Binding affinity value, ca. -10.00 kcal/mol, a Ki value 40 nM and Ligand

efficiency ca -0.15 kcal/mol. These results are comparable to the drug crystallized in the above- mentioned protein, currently used against TBC.

AUTHOR CONTRIBUTION

Ivan Vito Ferrari (I.V.F) conceived, designed and wrote the paper, performed the calculations and analysed the data. Paolo Patrizio (P.P) overview the article and wrote the paper.

DECLARATION OF COMPETING INTEREST

The authors declare they have no potential conflicts of interest to disclose

REFERENCES

1. Lin W, Mandal S, Degen D, Liu Y, Ebright YW, et al. Structural basis of Mycobacterium tuberculosis transcription and transcription inhibition. *Molecular cell*. 2017;66(2):169-79.
2. WHO. Tuberculosis (TB). 2020.
3. Trott O, Olson AJ. AutoDock Vina: improving the speed and accuracy of docking with a new scoring function, efficient optimization, and multithreading. *J comp chem*. 2010;31(2):455-61.
4. Morris GM, Huey R, Lindstrom W, Sanner MF, Belew RK, et al. AutoDock4 and AutoDockTools4: Automated docking with selective receptor flexibility. *J comp chem*. 2009;30(16):2785-91.
5. Pettersen EF, Goddard TD, Huang CC, Couch GS, Greenblatt DM, et al. UCSF Chimera—a visualization system for exploratory research and analysis. *J comp chem*. 2004;25(13):1605-12.
6. Sargis D. AutoDockTools, Pmv, and Vision. MGLTools. 2020.
7. Guex N, Peitsch MC. SWISS-MODEL and the Swiss-Pdb Viewer: an environment for comparative protein modeling. *electrophoresis*. 1997;18(15):2714-23.
8. Halgren TA. Merck molecular force field. I. Basis, form, scope, parameterization, and performance of MMFF94. *J. Comput. Chem*. 1996;17(5-6):490-519.
9. Biovia, Dassault Systèmes. "BIOVIA workbook, release 2017; BIOVIA pipeline pilot, release 2017." San Diego: Dassault Systèmes (2020).